

TRIVOREX[®]



POLYVALENT NEUTRALIZING ABSORBENT
TO FACILITATE THE MANAGEMENT OF LIQUID CHEMICAL SPILLS



- Trivorex[®] neutralises both acid and basic products and indicates the end of the reaction of neutralization by a colour change
- Trivorex[®] is an active absorbent, for all kinds of liquid spills: aqueous, hydrocarbons, oils, solvents...
- Trivorex[®] solidifies to facilitate the evacuation of liquid spills



PREVOR

ANTICIPATE AND SAVE

Toxicology Laboratory & Chemical Risk Management

Packaging



• 800 g POWDER SHAKER FOR LABORATORIES AND SMALL SPILLS

1



• Sprinkle to cover the puddle

2




• Contain the spill with Trivorex® in order to stop its spreading


3



4



• Wait for the colour to turn back to yellow-green* and for the product to solidify (if necessary add a little of water in case of highly concentrated products)
* indicator of a neutral pH



5



• Remove with a vacuum cleaner or a shovel
• Evacuate into the appropriate disposal device, in accordance with disposal regulations

Why use Trivorex®?

UNIDENTIFIED SPILL?

Trivorex® absorbs all liquids: clean and secure without a second thought!

CORROSIVE LIQUID SPILL?

Trivorex® minimizes the corrosive risk because it neutralizes both acidic and alkaline chemicals

HANDLING CHEMICALS IN LABORATORIES?

The Trivorex® shaker is a simple and essential device for handling safety in laboratories to absorb, solidify and evacuate accidental spills

LEAKAGE OF INCOMPATIBLE CHEMICALS AT THE BOTTOM OF YOUR TRASH BINS?


Trivorex® covers the bottom of your trash bins in order to soak up and isolate the various chemical species and prevent hazardous reactions from occurring

HANDLING OF HYDROFLUORIC ACID (HF)?

Trivorex® reduces the corrosivity and prevents the emission of harmful vapours

HAZARDOUS HANDLING OF AGGRESSIVE NEUTRALIZING ABSORBENTS?

Trivorex® is non toxic, non allergenic, non irritating



A LIQUID SPILL IN A DIFFICULT TO REACH AREA?

Because of its powder form, Trivorex® can reach into corners and will solidify to allow easy removal

USE OF RETENTION TANKS?

Trivorex® absorbs leftover liquids and prevents incompatible products from mixing together

A CMR OR TOXIC CHEMICAL SPILL?

Trivorex® encapsulates liquids and reduces the contact with harmful chemicals

A SPILL OF LIQUID EMITTING HARMFUL VOLATILE VAPOURS?

Trivorex® contains the emission of vapours by encapsulating liquids in a molecular network

OVERFLOWING LIQUIDS, LEAKING MACHINES?


Use Trivorex® to soak up runoff liquids

MULTIPLE HAZARDOUS SPILL SITES?

Trivorex® is a safety device adapted to any situation



• 10 Kg BUCKET FOR WORKSHOPS, AND PRODUCTION OR STORAGE AREAS



• BULK (200 Kg CONTAINER) TO MEET GREATER NEEDS

1



• Contain the spill with Trivorex® in order to stop its spreading

2




• Sprinkle to cover the puddle


3



• Wait for the colour to turn back to yellow-green* and for the product to solidify (if necessary add a little of water in case of highly concentrated products)
* indicator of a neutral pH




4



• Remove with a vacuum cleaner or a shovel
• Evacuate into the appropriate disposal device, in accordance with disposal regulation

5

Trivorex[®], a safety device adapted to any situation



Trivorex[®] is used and approved by the **AP-HP** (*Association of Parisian Hospitals*), and **A.F.A.Q.A.P.** (*French Association of Quality Assurance in Anatomy and Cytological Pathology*) and has been recommended by **NASA** (USA) and **Procter & Gamble** (Germany).

Trivorex[®] is used in numerous fields of activity: surface treatment, solar energy, chemical industry, hazardous material transportation, micro-electronics, metallurgy, glass industry, hospitals, food industry, water treatment, fire fighters, universities...

Trivorex[®] is non toxic, non allergenic, non irritating and non explosive
Trivorex[®] must be used in accordance with safety practices
and intervention measures in effect

PREVOR would be pleased to organize chemical risk management training sessions for your staff